

# ATLANTIS™

## Implant compatibility chart



Manufacturer	Implant	Titanium	GoldHue	Zirconia, white	Zirconia, shaded	Screwdriver type
Astra Tech	<b>Internal</b>					
	OsseoSpeed™/OsseoSpeed™ TX 3.0 - Yellow	x	x			0.050" Hex screwdriver
	OsseoSpeed™/OsseoSpeed™ TX 3.5, 4.0 - Aqua	x	x	x	x	
	OsseoSpeed™/OsseoSpeed™ TX 4.5, 5.0 - Lilac	x	x	x	x	
OsseoSpeed™ 4.0 - Blue (Previous)	x	x				
BioHorizons	<b>Internal</b>					
	Internal 3.5, 4.5, 5.7	x	x	x	x	1.25 mm Hex screwdriver
	Tapered Internal 3.5, 4.5, 5.7	x	x	x	x	
	<b>External</b>					
	External Implant 3.5	x	x			0.050" Hex screwdriver
	External Implant 4.0, 5.0	x	x	x	x	
Biomet 3i	<b>Internal</b>					
	Certain™ 3.25, 4.0, 4/3, 5.0, 5/4, 6.0	x	x	x	x	0.048" Hex screwdriver
	Certain™ Prevail™ 3/4/3, 4/5/4, 5/6/5, 4/3, 5/4, 6/5	x	x	x	x	
	Certain™ XP 4/5, 5/6	x	x	x	x	
	<b>External</b>					
	MicroMiniplant™ 3.25	x	x			0.048" Hex screwdriver
	Miniplant™ 3.25	x	x	x	x	
	Osseotite™ 3.75, 4.0, 5.0, 6.0	x	x	x	x	
Osseotite™ XP 4/5, 5/6	x	x	x	x		
Dentsply Friadent	<b>Internal</b>					
	XiVE S 3.0	x	x			0.048" Hex screwdriver
	XiVE S 3.4, 3.8, 4.5, 5.5	x	x	x	x	
Frialit 3.4, 3.8, 4.5, 5.5	x	x	x	x		
Keystone Dental	<b>External</b>					
	Restore™ Small Diameter 3.3, 3.4	x	x			0.048" Hex screwdriver
	Restore™ Regular Diameter 3.75, 4.0	x	x	x	x	
Restore™ Wide Diameter 5.0, 5.5, 6.0	x	x				
Nobel Biocare	<b>Internal</b>					
	NobelActive™ NP 3.5	x	x	x	x	Nobel Biocare Unigrip
	NobelActive™ RP 4.3, 5.0	x	x	x	x	
	NobelReplace™ NP 3.5	x	x	x	x	
	NobelReplace™ RP 4.0	x	x	x	x	
	NobelReplace™ WP 5.0	x	x	x	x	
	NobelReplace™ 6.0	x	x	x	x	
	<b>External</b>					
	Brånemark System™ NP 3.3	x	x			Nobel Biocare Unigrip
	Brånemark System™ RP 3.75, 4.0	x	x	x	x	
	Brånemark System™ WP 5.0, 5.5	x	x	x	x	
	SteriOss Replace™ 3.1, 3.5, 4.3, 5.0, 6.0	x	x			0.050" Hex screwdriver
SteriOss HL™ 3.25, 3.8, 4.5, 5.0, 6.0	x	x				
Straumann	<b>Internal</b>					
	Bone Level 3.3 NC (Narrow CrossFit™)	x	x			Straumann screwdriver compatible
	Bone Level 4.1, 4.8 RC (Regular CrossFit™)	x	x	x	x	
	Standard/Standard Plus 4.8 RN (Regular Neck)	x	x			
	Standard/Standard Plus 4.8 WN (Wide Neck)	x	x			
	<b>External</b>					
Standard Plus 3.3 NN (Narrow Neck)	x	x			Straumann screwdriver compatible	
Sybron Implant Solutions	<b>External</b>					
	Endopore 3.5, 4.1, 5.0 - external	x	x			0.050" Hex screwdriver
Entegra 3.5, 4.1, 5.0 - external	x	x				
Zimmer Dental	<b>Internal</b>					
	Tapered Screw-Vent™ 3.5, 4.5, 5.7	x	x	x	x	1.25 mm Hex screwdriver
	Screw-Vent™ 3.3, 3.7, 4.5	x	x	x	x	
SwissPlus™ 4.8	x	x				

## Turnaround times

Turnaround times (in-house business days) for all case sizes		
Titanium	Gold-shaded titanium	Zirconia
2	4	4

**Note:** Turnaround time listed begins the first business day after receipt of all case information and material required and does not include shipping, which is determined by your choice of delivery method. Turnaround for lab-based scanned cases is one day less than the times shown.

### Turnaround times are based on the following:

- Design approvals to process case must be granted by 10 am your local time.
- Lab-scanned cases, including rescans, must be uploaded to Atlantis™ WebOrder by 3 pm (EST) in order to save a day of processing compared with the times listed above.
- Material, on-hold and case information issues must be resolved by 8 pm (EST).

## Other services

Service	Additional turnaround days
Laboratory services	add 1 day
Provisional restoration	add 5 days
Removable diagnostic preparation	add 2 days

- For the most current implant compatibility, please visit [www.atlantisabutment.com](http://www.atlantisabutment.com) or contact Customer Service at 800-531-3481.
- All orders for Atlantis™ Abutments are shipped with the corresponding abutment screw needed at no extra charge.

**Note:** The implant manufacturer's recommended torque should always be used when placing the abutment screw.

