



G&H Dental Arts, Inc.

Connect Your CEREC® to G&H Dental Arts' Workflow



Save time and money when you send scans directly from your Cerec® digital impression scanner. **In addition to all-ceramic restorations, G&H accepts digital impressions for PFM and full-cast restorations.**

- Open and log into the Sirona Connect portal.
- Click "My Account" and then "My Favorite Laboratories".
- Search for G&H and click "Add" at the bottom of the page.

If you have any questions, contact G&H at 800-548-3384 or Cerec® at 800-944-6365 and ask for digital impression technical support.

Looking for more information on our products?

Save money on every digital submission by choosing model free when possible, visit gandhdental.com.

Digital model required on anterior and bridge cases at an additional fee.

SAVE \$50 On Your First Digital Case

Digital submissions are simple when you choose G&H Dental. Reference coupon code **GO_DIGITAL** on your Rx when you upload your first case and save \$50 on your first digital submission.

Coupon Code: GO_DIGITAL