



G&H Dental Arts, Inc.



Connect Your iTero® to G&H Dental Arts' Workflow

Save time and money when you send scans directly from your iTero® digital impression scanner. **In addition to all-ceramic restorations, G&H accepts digital impressions for PFM and full-cast restorations.**

- Call iTero® tech support at 800-577-8767, select option 2 and ask them to add G&H Dental Arts as a lab partner.
- Select G&H as your lab partner and submit your digital file.

If you have any questions, contact G&H at 800-548-3384 or iTero® at 800-577-8767 and ask for digital impression technical support.

Looking for more information on our products?

Save money on every digital submission by choosing model free when possible, visit gandhdental.com.

Digital model required on anterior and bridge cases at an additional fee.

SAVE \$50 On Your First Digital Case

Digital submissions are simple when you choose G&H Dental. Reference coupon code **GO_DIGITAL** on your Rx when you upload your first case and save \$50 on your first digital submission.

Coupon Code: GO_DIGITAL